


Anti-DlSt Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P36957

Item Number: CM0009383



Introduction

Anti-DLST rabbit monoclonal antibody for Western blot, IF/ICC, IHC-P, and ELISA. Validated in human samples. Detects ~49kDa protein, a component of the 2-oxoglutarate dehydrogenase complex involved in TCA cycle and histone succinylation.

[Learn More](#)

Parameter	Value
Product Name	DLST Rabbit mAb
Remarks/Alias	DLTS; KGD2; PGL7
Species	Human
GeneID (Human)	1743
GeneID	1743
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P36957
Protein Weight	49kDa
Shipping	Ice bag