

Anti-Dlx2 Rabbit Polyclonal Antibody For Wb, Elisa - Q07687

Item Number: CM0025720

Introduction

Anti-DLX2 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. DLX2 is a homeobox transcription factor involved in neuronal differentiation and craniofacial development.

[Learn More](#)

Parameter	Value
Product Name	DLX2 Rabbit pAb
Remarks/Alias	TES1; TES-1; DLX2
Species	Human
Gene ID	1746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human DLX2 (NP_004396.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q07687
Protein Weight	34kDa
Shipping	Ice bag