

Anti-Dlx3 (Ai4, Tdo) Rabbit Monoclonal Antibody For Wb, Elisa - O60479

Item Number: CM0002227



Introduction

Rabbit monoclonal antibody targeting human DLX3 (Homeobox protein DLX-3). Validated for WB and ELISA; cross-reacts with mouse. Ideal for transcription regulation and tooth development research.

[Learn More](#)

Parameter	Value
Product Name	DLX3 Rabbit mAb
Remarks/Alias	Ai4; TDO; DLX3
Species	Human
GeneID	1747
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 188-287 of human DLX3 (O60479)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60479
Protein Weight	32kDa
Shipping	Ice bag