

Anti-Dlx5 Polyclonal Antibody For Wb, Elisa - P56178

Item Number: CM0013656



Introduction

Rabbit polyclonal antibody against human DLX5, a BMP-responsive transcription factor critical for osteoblast differentiation and bone development. Suitable for WB and ELISA, cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	DLX5 Rabbit pAb
Remarks/Alias	SHFM1; SHFM1D; DLX5
Species	Human
GeneID (human)	1749
GeneID	1749
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human DLX5 (NP_005212.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56178
Protein Weight	32kDa
Shipping	Ice bag