

Anti-Dlx6 Polyclonal Antibody For Wb, Elisa - P56179

Item Number: CM0027384

Introduction

DLX6 rabbit polyclonal antibody validated for Western blot and ELISA applications. Detects human DLX6, a homeobox transcription factor involved in development. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	DLX6 Rabbit pAb
Remarks/Alias	DLX6
Species	Human
GeneID (human)	1750
GeneID	1750
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 119-293 of human DLX6 (NP_005213.3))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56179
Protein Weight	20kDa
Shipping	Ice bag