


Anti-Dmap1 Polyclonal Antibody For Wb, Elisa - Q9Npf5

Item Number: CM0019433



Introduction

Rabbit polyclonal antibody targeting DMAP1 (EAF2, SWC4), validated for Western blot and ELISA in human, mouse, rat. DMAP1 is a transcription regulator and chromatin modifier. 53 kDa.

[Learn More](#)

Parameter	Value
Product Name	DMAP1 Rabbit pAb
Remarks/Alias	EAF2; SWC4; MEAF2; DNMAP1; DNMTAP1; DMAP1
Species	Human
Gene ID (Human)	55929
Gene ID	55929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human DMAP1 (NP_061973.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPF5
Protein Weight	53kDa
Shipping	Ice bag