


Anti-Dmap1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Npf5

Item Number: CM0002176



Introduction

DMAP1 Rabbit Monoclonal Antibody suitable for WB and ELISA, targeting human DNA methyltransferase 1-associated protein 1 (Q9NPF5). Key component of the NuA4 histone acetyltransferase complex involved in transcription, DNA repair, and apoptosis. Bulk supply for IVD development.

[Learn More](#)

Parameter	Value
Product Name	DMAP1 Rabbit mAb
Remarks/Alias	EAF2; SWC4; MEAF2; DNMAP1; DNMTAP1; DMAP1
Species	Human
Gene ID (Human)	55929
Gene ID	55929
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NPF5
Protein Weight	53kDa
Shipping	Ice bag