


Anti-Dmc1 Polyclonal Antibody For Wb And Elisa - Q14565

Item Number: CM0026714



Introduction

Rabbit polyclonal antibody against human DMC1, a meiotic recombination protein. Validated for Western blotting and ELISA, cross-reacts with human, mouse, and rat. Suitable for meiosis and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	DMC1 Rabbit pAb
Remarks/Alias	DMC1H, LIM15, dj199H16.1, DMC1
Species	Human
Gene ID (Human)	11144
Gene ID	11144
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DMC1 (NP_008999.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14565
Protein Weight	38kDa
Shipping	Ice bag