


Anti-Dmd Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P11532

Item Number: CM0029745



Introduction

DMD Rabbit Polyclonal Antibody validated for Western blot (WB), immunofluorescence (IF-P), and ELISA in human, mouse, rat. Targets dystrophin, a 427 kDa cytoskeletal anchor essential for muscle sarcolemma integrity. Suitable for neuromuscular disease research.

[Learn More](#)

Parameter	Value
Product Name	DMD Rabbit pAb
Remarks/Alias	BMD; CMD3B; MRX85; DXS142; DXS164; DXS206; DXS230; DXS239; DXS268; DXS269; DXS270; DXS272
Species	Human
Gene ID (Human)	1756
Gene ID	1756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3414-3685 of human DMD (NP_003997.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11532
Protein Weight	427kDa
Shipping	Ice bag