

Anti-Dmp1 Polyclonal Antibody For Wb And Elisa - Q13316

Item Number: CM0021994

Introduction

CamelBio's anti-DMP1 polyclonal antibody (Q13316) is designed for WB and ELISA applications, targeting human, mouse, and rat DMP1. Ideal for osteoblast differentiation and biomineralization research.

[Learn More](#)

Parameter	Value
Product Name	DMP1 Rabbit pAb
Remarks/Alias	ARHP; ARHR; DMP-1; DMP1
Species	Human
GeneID (Human)	1758
GeneID	1758
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DMP1 (NP_004398.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13316
Protein Weight	56kDa
Shipping	Ice bag