

Anti-Dmp1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13316

Item Number: CM0016047

Introduction

Rabbit polyclonal antibody targeting human dentin matrix acidic phosphoprotein 1 (DMP1; also known as DMP-1, ARHP). Validated for WB and ELISA, cross-reacts with mouse. Ideal for osteoblast differentiation and biomineralization research.

[Learn More](#)

Parameter	Value
Product Name	DMP1 Rabbit pAb
Remarks/Alias	ARHP; ARHR; DMP-1; DMP1
Species	Human
Gene ID (Human)	1758
Gene ID	1758
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-513 of human DMP1 (NP_004398.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13316
Protein Weight	56kDa
Shipping	Ice bag