

Anti-Dmrt1 Polyclonal Antibody For Wb, Elisa - Q9Y5R6

Item Number: CM0025327



Introduction

Rabbit polyclonal antibody against human DMRT1 (39kDa), validated for WB and ELISA. Ideal for studying male sex determination, spermatogenesis, and testis biology.

[Learn More](#)

Parameter	Value
Product Name	DMRT1 Rabbit pAb
Remarks/Alias	DMT1; CT154; DMRT1
Species	Human
Gene ID (Human)	1761
Gene ID	1761
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 134-373 of human DMRT1 (NP_068770.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y5R6
Protein Weight	39kDa
Shipping	Ice bag