


Anti-Dmt1/Slc11A2 Monoclonal Antibody For Wb - P49281

Item Number: CM0004956



Introduction

Rabbit monoclonal antibody against DMT1/SLC11A2 (DMT1, NRAMP2), validated for WB, IF/ICC, IP, and ELISA. Reacts with human and mouse. Targets the divalent metal transporter essential for iron homeostasis.

[Learn More](#)

Parameter	Value
Product Name	DMT1/SLC11A2 Rabbit mAb
Remarks/Alias	DCT1; DMT1; AHMIO1; NRAMP2; DMT1/SLC11A2
Species	Human
Gene ID (Human)	4891
Gene ID	4891
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human DMT1/SLC11A2 (NP_000608.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49281
Protein Weight	61KDa/62KDa/64KDa/65KDa
Shipping	Ice bag