


Anti-Dmt1/Slc11A2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49281

Item Number: CM0017346



Introduction

Rabbit polyclonal antibody targeting human DMT1/SLC11A2 (NRAMP2), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for iron transport and homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	DMT1/SLC11A2 Rabbit pAb
Remarks/Alias	DCT1; DMT1; AHM1O1; NRAMP2; DMT1/SLC11A2
Species	Human
Gene ID (Human)	4891
Gene ID	4891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human DMT1/SLC11A2 (NP_001167597.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49281
Protein Weight	62kDa
Shipping	Ice bag