

Anti-Dna Ligase I Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18858

Item Number: CM0005218



Introduction

Rabbit monoclonal antibody against human DNA Ligase I, validated for Western blotting, immunohistochemistry, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Suitable for DNA repair, replication, and recombination research.

[Learn More](#)

Parameter	Value
Product Name	DNA Ligase I Rabbit mAb
Remarks/Alias	LIgI; IMD96; hLig1; DNA Ligase I
Species	Human
GeneID (Human)	3978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human DNA Ligase I (P18858).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18858
Protein Weight	102kDa
Shipping	Ice bag