

Anti-Dna Ligase Iii/Lig3 Monoclonal Antibody For Wb, Ihc-P, Elisa - P49916

Item Number: CM0006790



Introduction

Rabbit monoclonal antibody targeting human DNA Ligase III/LIG3 (P49916). Confirmed for WB, IHC-P, ELISA with human, mouse, rat reactivity. Ideal for DNA repair and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	DNA Ligase III/LIG3 Rabbit mAb
Remarks/Alias	LIG2; MTDPS20; LIG3alpha; DNA Ligase III/LIG3
Species	Human
GeneID (Human)	3980
GeneID	3980
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human DNA Ligase III/LIG3 (NP_039269.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49916
Protein Weight	113kDa
Shipping	Ice bag