

Anti-Dna Ligase Iii/Lig3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49916

Item Number: CM0013372



Introduction

Rabbit polyclonal antibody targeting DNA Ligase III (LIG3/LIG2) for WB, IHC-P, IF/ICC, and ELISA applications. Detects human, mouse, and rat proteins. Suitable for DNA repair, mitochondrial base-excision repair, and MTDP520 research.

[Learn More](#)

Parameter	Value
Product Name	DNA Ligase III/LIG3 Rabbit pAb
Remarks/Alias	LIG2; MTDP520; LIG3alpha; DNA Ligase III/LIG3
Species	Human
GeneID (Human)	3980
GeneID	3980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-460 of human DNA Ligase III/LIG3 (NP_039269.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49916
Protein Weight	113kDa
Shipping	Ice bag