

Anti-Dna Ligase Iii/Lig3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P49916

Item Number: CM0019806



Introduction

Rabbit polyclonal antibody targeting human DNA ligase 3 (LIG3) for WB, IP, ELISA. Reacts with human and mouse. Ideal for DNA repair and mitochondrial base-excision repair research. Molecular weight: 113 kDa.

[Learn More](#)

Parameter	Value
Product Name	DNA Ligase III/LIG3 Rabbit pAb
Remarks/Alias	LIG2; MTDPS20; LIG3alpha; DNA Ligase III/LIG3
Species	Human
GeneID (Human)	3980
GeneID	3980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 730-1009 of human DNA Ligase III/LIG3 (NP_039269.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49916
Protein Weight	113 kDa
Shipping	Ice bag