


Anti-Dna Ligase Iv Rabbit Monoclonal Antibody For Wb, Elisa - P49917

Item Number: CM0003250



Introduction

Rabbit monoclonal antibody targeting human DNA ligase IV (LIG4), validated for Western blot and ELISA. Essential for NHEJ DNA repair and V(D)J recombination studies, with a molecular weight of 104 kDa.

[Learn More](#)

Parameter	Value
Product Name	DNA Ligase IV Rabbit mAb
Remarks/Alias	LIG4S; DNA Ligase IV
Species	Human
Gene ID (Human)	3981
Gene ID	3981
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 504-685 of human DNA Ligase IV (P49917)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49917
Protein Weight	104kDa
Shipping	Ice bag