


Anti-Dna Ligase Iv Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49917

Item Number: CM0020581



Introduction

High-quality rabbit polyclonal antibody targeting human DNA Ligase IV (LIG4). Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for NHEJ and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	DNA Ligase IV Rabbit pAb
Remarks/Alias	LIG4S; DNA Ligase IV
Species	Human
Gene ID (Human)	3981
Gene ID	3981
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human DNA Ligase IV (NP_002303.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49917
Protein Weight	104kDa
Shipping	Ice bag