


Anti-Dna-Pkcs Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P78527

Item Number: CM0003774



Introduction

High-quality rabbit monoclonal antibody targeting human DNA-PKcs (SWISS: P78527). Validated for WB, IP, IF/ICC, IHC-P, and ELISA. Ideal for DNA damage repair and kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	DNA-PKcs Rabbit mAb
Remarks/Alias	HYRC; p350; DNAPK; DNPk1; HYRC1; IMD26; XRCC7; DNAPKc; DNA-PKC; DNA-PKcs
Species	Human
GeneID (Human)	5591
GeneID	5591
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 3993-4128 of human DNA-PKcs (NP_008835.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P78527
Protein Weight	469 kDa
Shipping	Ice bag