


Anti-Dna-Pkcs Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P78527

Item Number: CM0030923



Introduction

High-quality rabbit polyclonal antibody targeting human DNA-PKcs (PRKDC) protein. Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human and mouse. Ideal for DNA damage repair and kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	DNA-PKcs Rabbit pAb
Remarks/Alias	HYRC; p350; DNAPK; DNPK1; HYRC1; IMD26; XRCC7; DNAPKc; DNA-PKC; DNA-PKcs
Species	Human
GeneID (Human)	5591
GeneID	5591
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78527
Protein Weight	469kDa
Shipping	Ice bag