

Anti-Dna Polymerase Beta Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P06746

Item Number: CM0007541



Introduction

High-quality rabbit monoclonal antibody targeting DNA Polymerase beta, suitable for WB, IHC-P, and ELISA applications. Detects the 38kDa protein in human, mouse, and rat samples. A key reagent for base-excision repair research.

[Learn More](#)

Parameter	Value
Product Name	DNA Polymerase beta Rabbit mAb
Remarks/Alias	POLB; DNA polymerase beta; DNA Polymerase beta
Species	Human
GeneID (Human)	5423
GeneID	5423
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 236-335 of human DNA Polymerase beta (P06746)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06746
Protein Weight	38kDa
Shipping	Ice bag