

Anti-Dna Polymerase Beta Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06746

Item Number: CM0020779



Introduction

KO-validated rabbit polyclonal antibody targeting human DNA polymerase beta (POLB). Suitable for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse samples. Recombinant immunogen. Ideal for base-excision repair and DNA damage studies.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] DNA Polymerase beta Rabbit pAb
Remarks/Alias	POLB; ta
Species	Human
Gene ID (Human)	5423
Gene ID	5423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-335 of human DNA Polymerase beta (NP_002681.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06746
Protein Weight	38kDa
Shipping	Ice bag