


Anti-Dna Polymerase Delta Subunit 4 (Pold4/P12) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcu8

Item Number: CM0026168



Introduction

Rabbit polyclonal antibody targeting human DNA polymerase delta subunit 4 (POLD4/p12), suitable for Western blot (WB) and ELISA. Useful for investigating DNA replication, repair, and genomic stability.

[Learn More](#)

Parameter	Details
Product Name	POLD4 Rabbit pAb
Remarks/Alias	p12; POLD5; POLD4
Species	Human
Gene ID (Human)	57804
Gene ID	57804
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human POLD4 (NP_066996.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HCU8
Protein Weight	12kDa
Shipping	Ice bag