


Anti-Dna Polymerase Eta Polyclonal Antibody For Wb - Q9Y253

Item Number: CM0019777



Introduction

Rabbit polyclonal antibody against human DNA polymerase eta (POLH), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for translesion synthesis research.

[Learn More](#)

Parameter	Value
Product Name	DNA polymerase eta Rabbit pAb
Remarks/Alias	XPV; XP-V; RAD30; RAD30A; DNA polymerase eta
Species	Human
GeneID (Human)	5429
GeneID	5429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 322-520 of human DNA polymerase eta (NP_006493.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y253
Protein Weight	78kDa
Shipping	Ice bag