

Anti-Dna Polymerase Gamma Monoclonal Antibody For Wb, Elisa - P54098

Item Number: CM0002489

Introduction

Recombinant rabbit monoclonal antibody against human DNA Polymerase gamma (POLG). Validated for Western blot and ELISA, cross-reacts with mouse. Essential reagent for studies on mitochondrial DNA replication and repair.

[Learn More](#)

Parameter	Value
Product Name	DNA Polymerase gamma Rabbit mAb
Remarks/Alias	PEO; MDP1; SCAE; MIRAS; POLG1; POLGA; SANDO; MTDPS4A; MTDPS4B; DNA Polymerase gamma
Species	Human
Gene ID (Human)	5428
Gene ID	5428
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54098
Protein Weight	140 kDa
Shipping	Ice bag