


Anti-Dna Polymerase Gamma Rabbit Polyclonal Antibody For Wb, Elisa - P54098

Item Number: CM0027605



Introduction

High-quality rabbit polyclonal antibody against human DNA Polymerase gamma (POLG). Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for mitochondrial DNA replication and disease research.

[Learn More](#)

Parameter	Value
Product Name	DNA Polymerase gamma Rabbit pAb
Remarks/Alias	PEO; MDP1; SCAE; MIRAS; POLG1; POLGA; SANDO; MTDP54A; MTDP54B; DNA Polymerase gamma
Species	Human
GeneID (Human)	5428
GeneID	5428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-720 of human DNA Polymerase gamma (NP_002684.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54098
Protein Weight	140kDa
Shipping	Ice bag