

Anti-Dna Polymerase Kappa (Polk) Polyclonal Antibody For Wb And Elisa - Q9Ubt6

Item Number: CM0018713



Introduction

Rabbit polyclonal antibody against human DNA polymerase kappa (POLK, DINP, DINB1). Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Targets the translesion synthesis DNA repair protein.

[Learn More](#)

Parameter	Value
Product Name	POLK Rabbit pAb
Remarks/Alias	DINP; POLQ; DINB1; POLK
Species	Human
GeneID (Human)	51426
GeneID	51426
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human POLK (NP_057302.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBT6
Protein Weight	99kDa
Shipping	Ice bag