

Anti-Dna Topoisomerase I (Top1) Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P11387

Item Number: CM0015030



Introduction

Rabbit polyclonal antibody targeting human DNA topoisomerase I (TOP1) for WB, IHC-P, IP, and ELISA applications. Cross-reacts with mouse. Detects a 91 kDa protein involved in DNA topology and transcription.

[Learn More](#)

Parameter	Value
Product Name	DNA topoisomerase I (TOP1) Rabbit pAb
Remarks/Alias	TOPI; DNA topoisomerase I (TOP1)
Species	Human
GeneID (Human)	7150
GeneID	7150
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human TOP1 (NP_003277.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11387
Protein Weight	91kDa
Shipping	Ice bag