


Anti-Dna Topoisomerase II Alpha (Top2A) Polyclonal Antibody For If/Icc, Elisa - P11388

Item Number: CM0013197



Introduction

Rabbit polyclonal antibody targeting human DNA topoisomerase II alpha (TOP2A), validated for IF/ICC and ELISA. Suitable for cancer and cell biology research, recognizing a 174kDa protein. - P11388

[Learn More](#)

Parameter	Value
Product Name	DNA topoisomerase II alpha (TOP2A) Rabbit pAb
Remarks/Alias	TOP2; TP2A; TOP2A; TOP2alpha; DNA topoisomerase II alpha (TOP2A)
Species	Human
Gene ID	7153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1314-1503 of human DNA topoisomerase II alpha (TOP2A) (NP_001058.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P11388
Protein Weight	174kDa
Shipping	Ice bag