


Anti-Dna Topoisomerase II Alpha (Top2A) Rabbit Monoclonal Antibody For Wb - P11388

Item Number: CM0012077



Introduction

Rabbit monoclonal antibody against human DNA topoisomerase II alpha (TOP2A), validated for WB, IF/ICC, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for cell cycle, chromatin, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	DNA topoisomerase II alpha (TOP2A) Rabbit mAb
Remarks/Alias	TOP2; TP2A; TOPIIA; TOP2alpha; DNA topoisomerase II alpha (TOP2A)
Species	Human
GeneID (Human)	7153
GeneID	7153
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1432-1531 of human DNA topoisomerase II alpha (TOP2A) (P11388).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11388
Protein Weight	174kDa/178kDa/179kDa/183kDa
Shipping	Ice bag