

Anti-Dnaja1 Rabbit Polyclonal Antibody [Ko Validated] For Wb, If/Icc, Elisa - P31689

Item Number: CM0016910



Introduction

Rabbit polyclonal antibody against DNAJA1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for apoptosis, chaperone, and protein transport studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DNAJA1 Rabbit pAb
Remarks/Alias	DJ-2; DJA1; HDJ2; HSDJ; HSJ2; HSJ-2; HSPF4; NEDD7; hDJ-2; A1
Species	Human
GeneID	3301
Immunogen	Recombinant protein corresponding to amino acids 178-397 of human DNAJA1 (NP_001530.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31689
Protein Weight	45kDa
Shipping	Ice bag