


Anti-Dnaja3 Polyclonal Antibody For Wb, Elisa - Q96Ey1

Item Number: CM0029837



Introduction

Rabbit polyclonal antibody against DNAJA3 (TID1, hTid-1), validated for WB and ELISA. Recognizes human, mouse, rat DNAJA3. Ideal for apoptosis and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJA3 Rabbit pAb
Remarks/Alias	TID1; HCA57; hTID-1; DNAJA3
Species	Human
Gene ID (Human)	9093
Gene ID	9093
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 221-480 of human DNAJA3 (NP_005138.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EY1
Protein Weight	52kDa
Shipping	Ice bag