

## Anti-Dnajb11 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubs4

Item Number: CM0014266

### Introduction

Rabbit polyclonal antibody targeting human DNAJB11 (ERdj3), an ER co-chaperone. Validated for WB, IHC-P, ELISA, and cross-reacts with mouse and rat. Suitable for protein folding and ER stress studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | DNAJB11 Rabbit pAb                                                                                                  |
| Remarks/Alias    | DJ9; EDJ; Dj-9; ERJ3; PKD6; ABBP2; ERdj3; ERJ3p; ABBP-2; UNQ537; PRO1080; DNAJB11                                   |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 51726                                                                                                               |
| Gene ID          | 51726                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 89-358 of human DNAJB11 (NP_057390.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q9UBS4                                                                                                              |
| Protein Weight   | 41kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |