


Anti-Dnajb12 Rabbit Polyclonal Antibody For Wb & Elisa - Q9Nxx2

Item Number: CM0014141



Introduction

Rabbit polyclonal antibody against DNAJB12, a co-chaperone involved in protein folding, ERAD, and potassium channel maturation. Validated for WB and ELISA in human, mouse, and rat samples. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	DNAJB12 Rabbit pAb
Remarks/Alias	DJ10; DNAJB12
Species	Human
GeneID	54788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-409 of human DNAJB12 (NP_060096.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXW2
Protein Weight	42kDa
Shipping	Ice bag