

Anti-Dnajb14 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tbm8

Item Number: CM0027345

Introduction

DNAJB14 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human DNAJB14, a co-chaperone for protein quality control and ion channel maturation. Cross-reacts with mouse and rat. Ideal for ER stress and viral infection studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJB14 Rabbit pAb
Remarks/Alias	EGNR9427; PRO34683; DNAJB14
Species	Human
GeneID (Human)	79982
GeneID	79982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human DNAJB14 (NP_001026893.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBM8
Protein Weight	43kDa
Shipping	Ice bag