

Anti-Dnajb1/Hsp40 Monoclonal Antibody () For Intracellular Flow Cytometry - P25685

Item Number: CM0008471



Introduction

-conjugated rabbit monoclonal antibody targeting human DNAJB1/HSP40. Validated for intracellular flow cytometry, ideal for studying chaperone-mediated protein folding and heat shock response. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human DNAJB1/HSP40 mAb
Remarks/Alias	Hdj1; Sis1; HSPF1; Hsp40; RSPH16B
Species	Human
Gene ID (Human)	3337
Gene ID	3337
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 241-340 of human DNAJB1/HSP40 (NP_006136.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P25685
Protein Weight	27kDa/38kDa
Shipping	Ice bag