

Anti-Dnajb6 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - 075190

Item Number: CM0027078



Introduction

Rabbit polyclonal antibody targeting DNAJB6, a chaperone implicated in neurodegeneration and muscle disease. Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for protein aggregation studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJB6 Rabbit pAb
Remarks/Alias	DJ4, MRJ, DnaJ, HSJ2, HHDJ1, HSJ-2, MSJ-1, LGMD1D, LGMD1E, LGMDD1, DNAJB6
Species	Human
Gene ID (Human)	10049
Gene ID	10049
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 177-326 of human DNAJB6 (NP_490647.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075190
Protein Weight	36kDa
Shipping	Ice bag