


Anti-Dnajb9 Polyclonal Antibody For Wb And Elisa - Q9Ubs3

Item Number: CM0027557



Introduction

Rabbit polyclonal antibody against human DNAJB9 (ERdj4), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for endoplasmic reticulum chaperone and UPR research.

[Learn More](#)

Parameter	Value
Product Name	DNAJB9 Rabbit pAb
Remarks/Alias	MDG1, ERdj4, MDG-1, MST049, MSTP049, DNAJB9
Species	Human
Gene ID (Human)	4189
Gene ID	4189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-223 of human DNAJB9 (NP_036460.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBS3
Protein Weight	26kDa
Shipping	Ice bag