

Anti-Dnajc10 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ixb1

Item Number: CM0012965

Introduction

Rabbit polyclonal antibody against human DNAJC10 (ERdj5), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for studying ER protein quality control, ERAD, and stress signaling.

[Learn More](#)

Parameter	Value
Product Name	DNAJC10 Rabbit pAb
Remarks/Alias	JPDI; MTHr; ERdj5; PDIA19; DNAJC10
Species	Human
GeneID (Human)	54431
GeneID	54431
Immunogen	Recombinant fusion protein corresponding to amino acids 614-793 of human DNAJC10 (NP_061854.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IXB1
Protein Weight	91kDa
Shipping	Ice bag