

Anti-Dnajc12 Monoclonal Antibody For Wb, Elisa - Q9Ukb3

Item Number: CM0002756



Introduction

Rabbit monoclonal antibody against human DNAJC12, suitable for WB and ELISA; cross-reacts with mouse. A tool for investigating co-chaperone roles and hyperphenylalaninemia.

[Learn More](#)

Parameter	Value
Product Name	DNAJC12 Rabbit mAb
Remarks/Alias	JDP1; HPANBH4; DNAJC12
Species	Human
GeneID (Human)	56521
GeneID	56521
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 99-198 of human DNAJC12 (Q9UKB3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UKB3
Protein Weight	23kDa
Shipping	Ice bag