

Anti-Dnajc15 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Y5T4

Item Number: CM0012394

Introduction

Rabbit monoclonal antibody targeting human DNAJC15 (MCJ) for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for mitochondrial respiratory chain and protein transport studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJC15 Rabbit mAb
Remarks/Alias	MCJ; HSD18; DNAJD1; DNAJC15
Species	Human
GeneID (Human)	29103
GeneID	29103
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-150 of human DNAJC15 (Q9Y5T4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5T4
Protein Weight	16kDa
Shipping	Ice bag