

Anti-Dnajc18 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H819

Item Number: CM0021444

Introduction

Rabbit polyclonal antibody targeting human DNAJC18 (DnaJ homolog subfamily C member 18). Validated for WB and ELISA, recognizes human and rat protein. Involved in polyomavirus ER membrane penetration during infection.

[Learn More](#)

Parameter	Value
Product Name	DNAJC18 Rabbit pAb
Remarks/Alias	DNAJC18
Species	Human
Gene ID (Human)	202052
Gene ID	202052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-350 of human DNAJC18 (NP_689899.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H819
Protein Weight	42kDa
Shipping	Ice bag