

Anti-Dnajc19 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Da6

Item Number: CM0025825

Introduction

Rabbit polyclonal antibody against human DNAJC19, validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat DNAJC19 (TIM14). Suitable for mitochondrial import and cardiolipin remodeling studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJC19 Rabbit pAb
Remarks/Alias	PAM18; TIM14; TIMM14; DNAJC19
Species	Human
Gene ID (Human)	131118
Gene ID	131118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-116 of human DNAJC19 (NP_660304.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DA6
Protein Weight	12kDa
Shipping	Ice bag