


Anti-Dnajc2 Polyclonal Antibody For Wb, Elisa - Q99543

Item Number: CM0028636



Introduction

Rabbit polyclonal antibody against human DNAJC2 (ZRF1, MPP11). Validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for chaperone and chromatin regulation studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJC2 Rabbit pAb
Remarks/Alias	ZRF1; ZUO1; MPP11; MPHOSPH11; DNAJC2
Species	Human
GeneID (Human)	27000
GeneID	27000
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human DNAJC2 (NP_055192.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99543
Protein Weight	72kDa
Shipping	Ice bag