


Anti-Dnajc5 Polyclonal Antibody For Wb And Elisa - Q9H3Z4

Item Number: CM0018359



Introduction

Rabbit polyclonal antibody against human DNAJC5 (CSP/CLN4) for WB and ELISA. Cross-reacts with mouse and rat. Immunogen: recombinant protein (aa 1-198). For neuroscience and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	DNAJC5 Rabbit pAb
Remarks/Alias	CSP; NCL; CLN4; CLN4B; DNAJC5A; mir-941-2; mir-941-3; mir-941-4; mir-941-5; DNAJC5
Species	Human
Gene ID (Human)	80331
Gene ID	80331
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DNAJC5 (NP_079495.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H3Z4
Protein Weight	22kDa
Shipping	Ice bag