

Anti-Dnajc7 Polyclonal Antibody For Wb, Elisa - Q99615

Item Number: CM0029414

Introduction

Rabbit polyclonal antibody targeting DNAJC7, validated for Western blot and ELISA in human samples. Ideal for co-chaperone and steroid receptor folding research.

[Learn More](#)

Parameter	Value
Product Name	DNAJC7 Rabbit pAb
Remarks/Alias	DJ11; DJC7; TPR2; TTC2; DNAJC7
Species	Human
GeneID (Human)	7266
GeneID	7266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-380 of human DNAJC7 (NP_003306.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99615
Protein Weight	56kDa
Shipping	Ice bag