

Anti-Dnajc9 Polyclonal Antibody For Wb, Elisa - Q8Wxx5

Item Number: CM0020438



Introduction

Validate with Western blot and ELISA, this rabbit polyclonal antibody targets human DNAJC9, a dual histone and heat shock co-chaperone. Cross-reacts with mouse and rat; ideal for nucleosome assembly and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DNAJC9 Rabbit pAb
Remarks/Alias	JDD1; SB73; HDJC9; C9
Species	Human
GeneID (Human)	23234
GeneID	23234
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DNAJC9 (NP_056005.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXX5
Protein Weight	30kDa
Shipping	Ice bag