

Anti-Dnajc9 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wxx5

Item Number: CM0013943



Introduction

Anti-DNAJC9 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Targets human histone chaperone and HSP70 co-chaperone, ideal for nucleosome assembly and heat shock research.

[Learn More](#)

Parameter	Value
Product Name	DNAJC9 Rabbit pAb
Remarks/Alias	JDD1; SB73; HDJC9; DNAJC9
Species	Human
Gene ID (Human)	23234
Gene ID	23234
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DNAJC9 (NP_056005.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8WXX5
Protein Weight	30kDa
Shipping	Ice bag